INSERT REMOVAL & REPLACEMENT






- Leave model attached to bone clamp.
- Unscrew tibial shaft from tibial plug.
- 3 Release straps and reattach tibial shaft to tibial plug.
- Remove tibial plug from soft tissue cavity using the reattached tibial shaft as a handle.
- Semove meniscal cap from tibial plug.



Meniscal replacement:

- Orient the meniscus so the anterior portion matches that of the tibial plug.
- Install new meniscus by pressing it over tibial plug.
- 3 Reinstall plug and shaft in reverse order of removal.

CARE & MAINTENANCE



SAWBONES®

Should you have any questions, please contact:

SAWBONES


10221 SW 188th Street Vashon, WA 98070 (206) 463-5551 info@sawbones.com

SAWBONES EUROPE AB

Krossverksgatan 3 216 16 Malmö Sweden +46 40 650 70 00 info@sawbones.se

www.sawbones.com

Dry Knee Arthroscopy Model Item # 1401



Thank you for purchasing the Sawbones arthroscopy knee (dry model). This arthroscopy model may be used for diagnostic and visualization skills including triangulation, and pathology identification.

Menisectomy, loose body retrieval and some suturing applications may also be performed. This model may also be manipulated in full range of motion including varus and valgus positions. With a little care and maintenance your new model should provide you with many years of service.

This instructional brochure will take you through:

- * Components and anatomical structures
- * Set up instructions
- * Insert removal and replacement instructions
- * Product features
- * Care and maintenance
- * Available replacement parts





Components:

- Soft tissue knee
- 2 Tibial plug
- 3 Tibial shaft
- 4 Normal patella
- **6** Normal meniscus
- 6 Bone clamp
- C-clamp

(Note: Five different meniscal inserts and two patella inserts are available. Included are the normal anatomy meniscus and patella.)



If more space is desired within the joint

retaining straps.

space, remove one or both of the tibial plug

The arthroscopy knee model can be moved in flexion, extention, varus and valgus.



Knee soft tissue replacement Item # 1400-5



Meniscus, pathology Item # 1403



Tibial plug Item # 1116-21



Meniscus, pathology Item # 1404



Tibial shaft Item # 1116-22



Meniscus, pathology Item # 1405



Patella, normal Item # 1119-4




Meniscus, pathology Item # 1406



Patella fronds Item # 1407



Bone clamp Item # 1600



Meniscus, normal Item # 1402



C-clamp Item # 1601

Anatomical structures:

- * Normal anatomy patella
- * Proximal tibia
- * Collateral ligaments
- * Anterior/posterior ligaments
- * Normal anatomy meniscus
- * Synovial capsule

* Distal femur

SET UP



- Attach the bone and C-clamp combination to the table as shown.
- Secure the femoral shaft in the bone clamp.



This model comes with eight standard pre-molded portals. Additional portals may be created with a #11 scalpel and a sharp trochar.