

# INJECTION KNEE WITH AUDIO FEEDBACK

#1517-6

SET-UP INSTRUCTIONS

Thank you for purchasing the Sawbones conductive knee injection model.

This injection trainer can be used to teach or demonstrate needle placement for aspiration and joint injection. With the correct needle insertion, an audible tone will sound indicating proper placement of the chosen area. Palpation of the major anatomical landmarks may also be performed. Additional conductive structures may be added upon request. With a little care and maintenance your new model should provide you with many years of service.


SAWBONES®

#1517-6 and #1517-6-11

Stand is ordered separately



#1517-6



#1517-6-11



## COMPONENTS

#1517-6 Knee Injection Model include:



- Soft tissue injection knee
- Audio feedback box
- Needle assembly with lead
- Dry erase marker



Note: Replacement parts available at [www.sawbones.com](http://www.sawbones.com)

## OPERATION AND ANATOMICAL FEATURES

This Knee injection Model contain these anatomical points:



Locate the appropriate spot for injection (i.e. joint space bursa). Insert needle into location. The audio feedback box will emit a tone when the tip of the needle has made contact with the correct area.


Use the dry erase marker to identify injection targets.



## SET-UP INSTRUCTIONS

 **SAWBONES®**

**SET-UP INSTRUCTIONS  
#1517-6**



**Step 1:** Screw purple EMG needle onto syringe. Note that this is a standard Leuer lock connection.

**Step 2:** Plug the black male end pin plug from the syringe into the female black cord end coming from the knee.

**Step 3:** Plug the gray phone cord end into the receiving jack on the emitter.

**Step 4:** The battery is stored outside of the emitter box to extend life. Remove back cover from emitter box and install battery.

## REPLACEMENT PARTS



#1485-106

Knee Skin Sleeve  
Replacement




#1620-13

Conductive Needle Set,  
2" Needle length



#1620

Conductive audio  
feedback emitter



#1620-14

Needle Assembly with  
Lead, 2" Needle length

Note: Replacement parts available at  
[www.sawbones.com](http://www.sawbones.com)

## CARE & MAINTENANCE

### IMPORTANT

To maintain audio feedback, do not remove plunger assembly or needle from the syringe. This model is **NOT** intended for injection of fluids.

Your model can be wiped down with any mild cleaner. Use only the dry erase marker on the model. Standard ink can not be removed. Store in a dry, dark place out of direct contact with sunlight.

