

WRIST ARTHROSCOPY MODEL FOR CARPAL TUNNEL REPAIR

#1415-8

SET-UP INSTRUCTIONS

Thank you for purchasing the Sawbones wrist arthroscopy model for carpal tunnel repair.



This arthroscopy model may be used for various open and closed carpal tunnel release techniques. Visualization of the major anatomical structures may also be performed either with an arthroscope or in an open environment. Part: #1415-8 includes base, #1415-3 without base. With a little care and maintenance your new model should provide you with many years of service.



#1415-8

COMPONENTS & ANATOMY

#1415-8 Wrist Arthroscopy Model for Carpal Tunnel repair include:

- Soft tissue hand
- Skin cover
- Transverse ligament
- Palmar insert
- Base
- Hyperextension foam cylinder



Anatomical structure:

- Carpal metacarpal bones
- Ulna
- Radius
- Flexor tendons
- Median nerve
- Ulnar artery
- Transverse ligament




#1415-3

#1415-3 Wrist Arthroscopy Model for Carpal Tunnel repair (without base) include:

- Soft tissue hand
- Skin cover
- Transverse ligament
- Palmar insert

Note: Replacement parts available at www.sawbones.com

REPLACEMENT PARTS



#1415-8-2

Transverse wrist ligament

#1415-8-1

Palmar insert

