

RediStik™ Wearable Task Trainer Products




Versatile, easy set-up simulation training kits that provide wearable, durable and realistic task trainers.

Simulation task trainers designed and vetted by frontline innovators at Texas Children's Hospital® in collaboration with Sawbones®.

RediStik™ Wearable PIV Trainer Kit

- **Self-sealing textured skin** and **subcutaneous tissue** provide the ability to practice:
 - Skin preparation techniques with alcohol or CHG application.
 - Visualization, palpation, and puncture of the vein with blood flash confirmation.
 - Threading of the catheter into the intravascular space.
 - Drawing back blood and administering saline flush.
 - Dressing placement.
- **Lightweight** and **versatile model** with a contoured design that quickly can be secured to a person or mannequin's arm as well as used in tabletop simulations.
- **Needle-proof backing** prevents accidental needle sticks.
- **One small and one large vein option** accommodates 18 GA to 24 GA needles.
- **Clinical testing demonstrated life expectancy averages 1,000+ punctures with 24 GA needles.**



Quick setup for just-in-time training and simulation scenarios.



Progress from novice to advanced in venipuncture and PIV insertion.


Order references:

- | | |
|----------------|---|
| ■ #1703-459 | RediStik™ Wearable PIV Trainer Kit, Light Skin Tone |
| ■ #1703-459-10 | RediStik™ Wearable PIV Trainer Kit, Dark Skin Tone |




RediStik™ Wearable Port Trainer Kit

- **Textured skin** and **subcutaneous tissue** provide the ability to practice:
 - Sterile skin preparation techniques with CHG application.
 - Life-like palpation and puncture of the titanium port.
 - Confirmation of needle placement with blood aspiration and ability to flush saline.
 - Dressing application.
- Includes simulation options for **superficial, normal, deep, angled, and float** port positioning depths.
- **Lightweight and versatile model** with a contoured design that quickly can be secured to a person or mannequin, with the attached shoulder harness system, as well as used in tabletop simulations.
- **Needle proof backing** prevents injury when worn on a person.
- Angiodynamics® Medium **Titanium Smart Port** option incl. in the kit.
- **Clinical testing demonstrated life expectancy for removable skin averages 1,200+ punctures.**



Worn on a person or mannequin to increase staff engagement.



Progress from novice to advanced in port access, de-access, positioning and trouble shooting.

Order references:

- | | |
|----------------|--|
| ■ #1703-501-10 | RediStik™ Wearable Port Trainer Kit with Titanium Smart Port |
| ■ #1703-501 | RediStik™ Wearable Port Trainer Kit without Port |





in collaboration with



RediStik™ Wearable CVC Trainer Kit

- **Wearable design, resilient skin, and pre-inserted Double Lumen CVC** provide the ability to practice:
 - Sterile technique, **skin preparation** with alcohol or chlorhexidine gluconate (CHG) and dressing application.
 - **Collecting blood cultures, lab draws, and flushing.**
 - Cap changes and trouble shooting.
 - Wearable design helps learners practice positioning techniques.
- Trainer provides opportunity to **practice institutional policy for CVC care.**
- **Novel material designed to provide repeated skin preparation and dressing applications. Dressing can remain on while in storage.** Dressings do not leave adhesive residue.
- **Lightweight and contoured design** allows for quick securement to a person's or mannequin's arm.
- **Versatile model** allows for either tabletop or wearable simulation options.



K
E
Y

F
E
A
T
U
R
E
S



Practice institutional policy for CVC care.

Order references:

- | | |
|----------------|---|
| ■ #1703-873 | RediStik™ Wearable CVC Trainer Kit, Light Skin Tone |
| ■ #1703-873-10 | RediStik™ Wearable CVC Trainer Kit, Dark Skin Tone |

To purchase RediStik™ products contact Sawbones Customer Service (contact details below).

Visit texaschildrens.org/redistik to access product details, skills video tips and techniques and data outcomes.



in collaboration with

